

CUSTOMER SERVICE BULLETIN

Radio Shack® CUSTOMER SERVICE BULLETIN

400 Atrium—One Tandy Center, Fort Worth, Texas 76102

Bulletin No.	1568-03
Description	Medical Office System
Model	Model III
Stock No.	26-1568 Version 1.2
Date	Dec. 22, 1981
Page 1 of	1

Problem

The date does not change from the initial date.

Solution

The problem is corrected by following the steps listed below.

1. Backup the diskette(s) and make the changes on the Backup copy of the program.
2. In BASIC load the program by typing LOAD "DLYJOUR".
3. Make the following corrections:

CHANGES (Retype the line or refer to the Edit section of the owners manual)

Old Line: 1160 LPRINT USING "TOTAL RECEIPTS: #####.##"; TT#(2)+TT#(3): LPRINT USING "TOTAL CHARGES: #####.##"; TT#(1): FORE=1 TO 3: TT#(E)=0: NEXTE: LPRINT CHR\$(12); : POKE 16425, 1: IFFQANDMH<=LH THEN D8\$=DH\$: CN\$="": GOSUB 1055 : GOTO 1080 ELSE IF NQ>0 AND FQ THEN 1200 ELSE 560

New Line: 1160 LPRINT USING "TOTAL RECEIPTS: #####.##"; TT#(2)+TT#(3): LPRINT USING "TOTAL CHARGES: #####.##"; TT#(1): FORE=1 TO 3: TT#(E)=0: NEXTE: LPRINT CHR\$(12); : POKE 16425, 1: IFFQANDMH<=LH THEN D8\$=DH\$: GOSUB 340: CN\$="": GOSUB 1055: GOTO 1080 ELSE IF NQ>0 AND FQ THEN 1200 ELSE 560

4. Type SAVE "DLYJOUR" to save the changes in the program.
5. At TRSDOS Ready, make a backup copy of the corrected diskette.